

KAT2A Rabbit pAb

货号: **B13702**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IF/ICC IHC
推荐浓度	WB: 1:100 - 1:500 IF/ICC: 1:50 - 1:200 IHC: 1:50 - 1:200
理论分子量	94kDa
实测分子量	94KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human KAT2A (NP_066564.2).
序列	MAEPSQAPTAPAAQPRPLQSPAPAPTPTPAPSPASAPIPTTPAPAPAPAAAPAGSTGTGGPGVGSGGAGSGGDPARPG LSQQQRASQRKAQVRGLPRA

靶点信息

研究背景	KAT2A, or GCN5, is a histone acetyltransferase (HAT) that functions primarily as a transcriptional activator. It also functions as a repressor of NF-kappa-B (see MIM 164011) by promoting ubiquitination of the NF-kappa-B subunit RELA (MIM 164014) in a HAT-independent manner (Mao et al., 2009 [PubMed 19339690]). [supplied by OMIM, Sep 2009]
------	--

基因ID	2648
基因名	KAT2A
Swiss	Q92830
别名	KAT2A;GCN5;GCN5L2;PCAF-b;hGCN5

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn