

CUX1 Rabbit pAb

货号: B13704

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus IHC: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	72-77kDa/147-165kDa
实测分子量	200KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,U-251MG,Jurkat,Mouse heart
细胞定位	Golgi apparatus membrane,Single-pass type IV membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CUX1 (NP_001189472.1).
序列	MAANVGSMFQYWKRFDLQQLQRELDATATVLANRQDESEQSRKRLIEQSREFKKNTPEDLRKQVAPLLKSFQGEIDALS KRSKEAAEAFNLNVYKRLIDVPDPVPALDLGQQLQLKVQRLHDIETENQKLRETLEEYNKEFAEVKNQEVTIKALKEKIREYEQ TLKNQAETIALEKEQKLQNDFAEKERKLQETQMSTTSKLEEAHKVQSLQTALEKTRTELFDLKTKYDEETTAKADEIEMIMT DLERANQRAEVAQREAETLREQLSSANHSLQLASQ

靶点信息

研究背景	The protein encoded by this gene is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	1523
基因名	CUX1
Swiss	P39880,Q13948
别名	CUX1;CASP;CDP;CDP/Cut;CDP1;COY1;CUTL1;CUX;Clox;Cux/CDP;GOLIM6;Nbla10317;p100;p110;p200;p75

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn