

usp29 Rabbit pAb

货号: B13739

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	104kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,SH-SY5Y,U-251MG
细胞定位	cytosol,nucleus,perinuclear region of cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human usp29 (NP_065954.1).
序列	LKYIQSNRKNPSSLEDLEKDRDLKLGPSFNTNCNGNPNLDETVLATQTLNAKNGLTSPLEPEHSQGDPRCNKAQVPLDSHSQQLQQGFPNLGNTCYMNAVL

靶点信息

研究背景	Predicted to enable cysteine-type endopeptidase activity and thiol-dependent deubiquitinase. Predicted to be involved in G1/S transition of mitotic cell cycle and protein deubiquitination. Predicted to be active in cytosol and nucleus.
基因ID	57663

基因名	USP29
Swiss	Q9HBJ7
别名	USP29;HOM-TES-84/86;usp29

产品验证



Western blot analysis of usp29 expressed in U-87MG,SH-SY5Y,U-251MG using usp29 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn