

CAPZA2 Rabbit pAb

货号: B13756

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Rattus norvegicus WB: Rattus norvegicus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	19kDa/32kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,THP-1,22Rv1,Mouse brain,Mouse heart
细胞定位	actin cytoskeleton,cortical cytoskeleton,cytosol,extracellular exosome,extracellular region
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-286 of human CAPZA2 (NP_006127.1).
序列	MADLEEQLSDEEKVRIAAKFIIHAPPGEFNEVFNDVRLLLNNDNLLREGAAHAFAYNLDQFTPVKIEGYEDQVLITEHGD L GNGKFLDPKNRICFKFDHLRKEATDPRPCEVENAVESWRTSVETALRAYVKEHYPNGVCTVYGKKIDGQQTIIACIESHQ F QAKNFWNGRWRSEWKFTITPSTTQVVGILKIQVHYEDGNVQLVSHKDIQDSLTVSNEVQTAKFEKIVEAAENEYQTAISE NYQTMSDTTFKALRRQLPVTRTKIDWNKILSYKIGKEMQNA

靶点信息

研究背景	The protein encoded by this gene is a member of the F-actin capping protein alpha subunit family. It is th e alpha subunit of the barbed-end actin binding protein Cap Z. By capping the barbed end of actin filame nts, Cap Z regulates the growth of the actin filaments at the barbed end.
------	---

基因ID	830
基因名	CAPZA2
Swiss	P47755
别名	CAPZA2;CAPPA2;CAPZ

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn