

# CSF2RA Rabbit pAb

货号: B13761

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                            |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  | <b>WB:</b> Mus musculus                                                                                |
| 应用    | <a href="#">WB</a>                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                              |
| 理论分子量 | 26kDa/31kDa/32kDa/38kDa/43kDa/46kDa/49kDa                                                              |
| 实测分子量 | 80kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | HeLa,DU145,HepG2,Mouse spleen,Mouse lung                                                               |
| 细胞定位  | Cell membrane,Secreted,Single-pass type I membrane protein                                             |
| 纯化    | Affinity purification                                                                                  |


抗原信息

|      |                                                                                                                                                                                                                                                                                                                        |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 23-320 of human CSF2RA (NP_758448.1).                                                                                                                                                                                                    |
| 序列   | EKSDLRTVAPASSLNVRFDSRTMNLSWDCQENTTFSKCFLTDKKNRVVEPRLSNNECSCTFREICLHEGVTFEVHVNTSQ<br>RGFQQKLLYPNSGREGTAAQNFSCFIYNADLMNCTWARGPTAPRDVQYFLYIRNSKRRREIRCPYYIQDSGTHVGCHLDN<br>LSGLTSRNYFLVNGTSREIGIQFFDSLDDTKKIERFNPPSNVTVCNTTHCLVRWKQPRTYQKLSYLDFQYQLDVHRKNTQ<br>PGTENLLINVSGDLENRYNFPSSSEPRAKHSVKIRAADVRLNWSSWSEAEIEFGSDDG |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | The protein encoded by this gene is the alpha subunit of the heterodimeric receptor for colony stimulating factor 2, a cytokine which controls the production, differentiation, and function of granulocytes and macrophages. The encoded protein is a member of the cytokine family of receptors. This gene is found in the pseudoautosomal region (PAR) of the X and Y chromosomes. Multiple transcript variants encoding different isoforms have been found for this gene, with some of the isoforms being membrane-bound and others being soluble. |
| 基因ID  | 1438                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 基因名   | CSF2RA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Swiss | P15509                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 别名    | CSF2RA;CD116;CDw116;CSF2R;CSF2RAX;CSF2RAY;CSF2RX;CSF2RY;GM-CSF-R-alpha;GMCSFR;GMR;SMDP4                                                                                                                                                                                                                                                                                                                                                                                                                                                                |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)