

SCNN1B Rabbit pAb

货号: B13792

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	72kDa/77kDa
实测分子量	75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,LO2,Mouse lung
细胞定位	Apical cell membrane,Cytoplasmic vesicle membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 406-640 of human SCN N1B (NP_000327.2).
序列	YPLPRGEKYCNNRDFPDWAHCYSDLQMSVAQRETCIGMCKESCNDTQYKMTISMADWPSEASEDWIFHVLSQERDQST NITLSRKGIVKLNIFYQEFNYRTIEESAANNIVWLLSNLGGQFGFWMGGSVLCLEFGEIIIDFVWITIILKVALAKSLRQRAQ ASYAGPPPTVAELVEAHTNFGFQPDTPRSPNTGPYPSEQALPIPGTPPPNYDSLRLQPLDVIESDSEGDAI

靶点信息

研究背景	Nonvoltage-gated, amiloride-sensitive, sodium channels control fluid and electrolyte transport across epit helia in many organs. These channels are heteromeric complexes consisting of 3 subunits: alpha, beta, a nd gamma. This gene encodes the beta subunit, and mutations in this gene have been associated with ps eudohypoaldosteronism type 1 (PHA1), and Liddle syndrome.
基因ID	6338

基因名	SCNN1B
Swiss	P51168
别名	SCNN1B;BESC1;ENaCb;ENaCbeta;SCNEB

产品验证



Western blot analysis of SCNN1B expressed in 293T,LO2,Mouse lung using SCNN1B Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. B locking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn