

NR0B1 Rabbit pAb

货号: B13793

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	CHIP-seq: Arabidopsis thaliana
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:10 - 1:100
理论分子量	43kDa/51kDa
实测分子量	48KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,NCI-H460,MCF7
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human NR0B1 (NP_000466.2).
序列	MAGENHQWQGSILYNMLMSAKQTRAAPEAPETRLVDQCWGCSCGDEPGVGREGLLGGRNVALLYRCCFCGKDHPRQG SILYSMLTSAKQTYAAPKAPEATLGPCWGCSCGSDPGVGRAGLPGGRPVALLYRCCFCGEDHPRQGSILYSLLTSSKQTH VAPAAPEARPGGAWWDRSYFAQRPGGKEALPGGRATALLYRCCFCGEDHPQQGSTLYCVPTSTNQAQAAPEERPRAPW WDTSSGALRPVALKSPQVVCEAASA

靶点信息

研究背景	This gene encodes a protein that contains a DNA-binding domain. The encoded protein acts as a dominant-negative regulator of transcription which is mediated by the retinoic acid receptor. This protein also functions as an anti-testis gene by acting antagonistically to Sry. Mutations in this gene result in both X-linked congenital adrenal hypoplasia and hypogonadotropic hypogonadism.
------	---

基因ID	190
基因名	NR0B1
Swiss	P51843
别名	NR0B1;AHC;AHCH;AHX;DAX-1;DAX1;DSS;GTD;HHG;NR0B1;SRXY2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn