

PTH1R Rabbit pAb

货号: B13794

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	66kDa
实测分子量	66kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,MCF7,U-251MG,DU145,THP-1,Mouse kidney,Mouse liver
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 464-593 of human PTH1 R (NP_000307.1).
序列	GEVQAEIKKSWSRWTLALDFKRKARSGSSSYSGPMVSHTSVTNVGPVGLGLPLSPRLPTATTNGHPQLPGHAKPGT PALETLETPPAMAAPKDDGFLNGSCSGLDEEASGPERPPALLQEEWETVM

靶点信息

研究背景	The protein encoded by this gene is a member of the G-protein coupled receptor family 2. This protein is a receptor for parathyroid hormone (PTH) and for parathyroid hormone-like hormone (PTHrP). The activity of this receptor is mediated by G proteins which activate adenylyl cyclase and also a phosphatidylinositol-calcium second messenger system. Defects in this receptor are known to be the cause of Jansen's metaphyseal chondrodysplasia (JMC), chondrodysplasia Blomstrand type (BOCD), as well as enchondromatosis. Two transcript variants encoding the same protein have been found for this gene.
------	--

基因ID	5745
基因名	PTH1R
Swiss	Q03431
别名	PTH1R;PFE;PTHR;PTHR1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn