

IFNAR1 Rabbit pAb

货号: B13798

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Mus musculus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	49kDa/56kDa/63kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,K562
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 470-557 of human IFNA R1 (NP_000620.2).
序列	PSLKPSSSIDEYFSEQPLKNLLSTSEEQIEKCFIENISTIATVEETNQTDDEHKKYSSQTSQDSGNYSNEDESESKTSEELQ QDFV

靶点信息

研究背景	The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinase s, which in turn phosphorylate several proteins, including STAT1 and STAT2. The encoded protein also functions as an antiviral factor.
基因ID	3454

基因名	IFNAR1
Swiss	P17181
别名	IFNAR1;AVP;IFN-alpha-REC;IFNAR;IFNBR;IFRC

产品验证



Western blot analysis of IFNAR1 expressed in SW480,K562 using IFNAR1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn