

GNAT3 Rabbit pAb

货号: B13837

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rattus norvegicus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	40kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-130 of human GNAT 3 (NP_001095856.1).
序列	GYSEQECMEFKAVIYSNTLQSI LAIVKAMTTLGIDYVNP R SAEDQRQLYAMANTLEDGGMTPQLAEVIKRL

靶点信息

研究背景	Sweet, bitter, and umami tastes are transmitted from taste receptors by a specific guanine nucleotide binding protein. The protein encoded by this gene is the alpha subunit of this heterotrimeric G protein, which is found not only in the oral epithelium but also in gut tissues. Variations in this gene have been linked to metabolic syndrome.
基因ID	346562

基因名	GNAT3
Swiss	A8MTJ3
别名	GNAT3;GDCA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn