

PLVAP Rabbit pAb

货号: B13844

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	50kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat lung,Rat kidney
细胞定位	Cell membrane,Cytoplasm,Membrane,Single-pass type II membrane protein,caveola,perinuclear region
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 63-442 of human PLVAP (NP_112600.1).
序列	TERRAEGLYSQLLGLTASQSNLTKELNFTTRAKDAIMQMWLNARRDLDRINASFRQCQGDRVIYTNNQRYMAAILSEKQC RDQFKDMNKSCDALLFMLNQKVKTLEVEIAKEKTICTDKESVLLNKRVAEEQLVECVKTRELQHQERQLAKEQLQKVQAL CLPLDKDKFEMDLRNLWRDSIIPRSLDNLGYNLYHPLGSELASIRRACDHMPSLMSSKVEELARSLRADIERVARENSDLQ RQKLEAQQGLRASQEAKQKVEKEAQAREAKLQAECSRQTQLALEEKAVLRKERDNLAKELEEKKREAEQLRMELAIRNSAL DTCIKTKSQPMMPVSRPMGPVNPQPIDPASLEEFKRKILESQRPPAGIPVAPSSG

靶点信息

研究背景	Endothelial cell-specific membrane protein involved in the formation of the diaphragms that bridge endothelial fenestrae. It is also required for the formation of stomata of caveolae and transendothelial channels . Functions in microvascular permeability, endothelial fenestrae contributing to the passage of water and solutes and regulating transcellular versus paracellular flow in different organs. Plays a specific role in embryonic development.
基因ID	83483
基因名	PLVAP
Swiss	Q9BX97
别名	PLVAP;FELS;PV-1;PV1;gp68

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn