

TEX11 Rabbit pAb

货号: B13851

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	71kDa/107kDa/108kDa
实测分子量	108kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	DU145,Mouse testis,Rat testis
细胞定位	Chromosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 701-940 of human TEX11 (Q8IYF3).
序列	AFEQTMFLSRALEEIQTCNDIHNFLKQTGTFSNDSCEKLLLLYEFVRAKLNDPLLESFLESVWELPHLETKTFETIAIIAMEKP AHYPLIALKALKKALLYKKEEPIIDISQYSKCMHNLVNLSVPDGASNVELCPLEEVWGYFEDALSHISRTKDYPEMILWLMV KSWNTGVLMFSRSKYASAEKWCGLALRFLNHLTSFKESYETQMNMLYSQLVEALSNNKGPVFHEHGYWSKSD

靶点信息

研究背景	This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variant s encoding distinct isoforms have been found for this gene.
基因ID	56159

基因名	TEX11
Swiss	Q8IYF3
别名	TEX11;SPGFX2;TGC1;TSGA3

产品验证



Western blot analysis of TEX11 expressed in DU145,Mouse testis,Rat testis using TEX11 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn