

ZNF423 Rabbit pAb

货号: B13855

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rattus norvegicus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	137kDa/144kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse eye,Rat eye
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 220-450 of human ZNF423 (NP_055884.2).
序列	KPFKCTVCKRGFSSTSSLQSHMQAHKKNEHLAKSEKEAKKDDFMCDYCEDTFSQTEELEKHVLTRHPQLSEKADLQCIHCPEVFDENTLLAIIHQAHANQKHKCPMCPEQFSSVEGVYCHLDSHRQPDSSNHSVSPDPVLGVSASMSSATPDSSASVERGSTPDSTLKPLRGQKKMRDDGQGWTKVVYSCPYCSKRDFNSLAVLEIHLKTIHADKPQQSHTCQICLDSM

靶点信息

研究背景	The protein encoded by this gene is a nuclear protein that belongs to the family of Kruppel-like C2H2 zinc finger proteins. It functions as a DNA-binding transcription factor by using distinct zinc fingers in different signaling pathways. Thus, it is thought that this gene may have multiple roles in signal transduction during development. Mutations in this gene are associated with nephronophthisis-14 and Joubert syndrome-19. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	23090
基因名	ZNF423
Swiss	Q2M1K9
别名	ZNF423;Ebfaz;JBTS19;NPHP14;OAZ;Roaz;ZFP423;Zfp104;hOAZ

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn