

KIF7 Rabbit pAb

货号: B13878

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	IP: Drosophila melanogaster WB: Drosophila melanogaster
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	150kDa
实测分子量	150kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,LO2,A-549
细胞定位	Cell projection,cilium
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1044-1343 of human KI F7 (NP_940927.2).
序列	QLDEAIEALDAAIEYKNEAITCRQRVLRASASLLSQCEMNLMAKLSYLSSETRALLCKYFDKVVTLREEQHQQQIAFSELEM QLEEQQRLVYWLEVALERQRLEMDRQLTLQQKEHEQNMQLLLQQSRDHLGEGLADSRQYEARIQALEKELGRYMWINQ ELKQKLGGVNAVGHSRGGEKRSKLCSEGRQAPGNEDELHLAPELLWLSPLTEGAPRTREETRDLVHAPLPLTWKRSSLCGE EQGSPEELRQREAAEPLVGRVLPVGEAGLPWNFGPLSKPRRELRRASPGMIDVRKNPL

靶点信息

研究背景	This gene encodes a cilia-associated protein belonging to the kinesin family. This protein plays a role in the sonic hedgehog (SHH) signaling pathway through the regulation of GLI transcription factors. It functions as a negative regulator of the SHH pathway by preventing inappropriate activation of GLI2 in the absence of ligand, and as a positive regulator by preventing the processing of GLI3 into its repressor form. Mutations in this gene have been associated with various ciliopathies.
基因ID	374654
基因名	KIF7
Swiss	Q2M1P5
别名	KIF7;ACLS;AGBK;HLS2;JBTS12;UNQ340

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn