

HOOK3 Rabbit pAb

货号: B13881

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	83kDa
实测分子量	170kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa,293T,mouse heart,mouse kidney,rat lung
细胞定位	Cytoplasm,Golgi apparatus,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 489-718 of human HOO K3 (NP_115786.1).
序列	IALLQSLLDDANLRKNELETENRLVNQRLLLEVQSQVEELQKSLQDQGSKAEDSVLLKKKLEEHLEKLHEANNELQKKRAIIE DLEPRFNSSLKIEELQEALRKKEEEMKQMEERYKKYLEKAKSVIRTLDPKQNQGAAPEIQALKNQLQERDRLFHSLEKEYE KTKSQREMEEKYIVSAWYNMGMTLHKKAAEDRLASTGSGQSFLARQRQATSSRRSYPGHVQPATAR

靶点信息

研究背景	Hook proteins are cytosolic coiled-coil proteins that contain conserved N-terminal domains, which attach to microtubules, and more divergent C-terminal domains, which mediate binding to organelles. The Drosophila Hook protein is a component of the endocytic compartment.
基因ID	84376

基因名	HOOK3
Swiss	Q86VS8
别名	HOOK3;HK3

产品验证



Western blot analysis of HOOK3 expressed in mouse heart, rat lung, LO2 using HOOK3 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn