

UNC93B1 Rabbit pAb

货号: B13899

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	RT-PCR: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	66kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,Mouse kidney,Mouse lung,Mouse spleen,Rat kidney,Rat lung,Rat spleen
细胞定位	Cytoplasmic vesicle,Endoplasmic reticulum membrane,Endosome,Lysosome,Multi-pass membrane protein,phagosome
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 450-550 of human UNC93B1 (NP_112192.2).
序列	KTGLSTLLGILYEDKERQDFIFTIYHWWQAVAIFTVYLGSSLHMKAKLAVLLVTLVAAAVSYLRMEQKLRRGVAPRQPRIPRPQHKVRGYRYLEEDNSDES

靶点信息

研究背景	This gene encodes a protein that is involved in innate and adaptive immune response by regulating toll-like receptor signaling. The encoded protein traffics nucleotide sensing toll-like receptors to the endolysosome from the endoplasmic reticulum. Deficiency of the encoded protein has been associated with herpes simplex encephalitis.
基因ID	81622

基因名	UNC93B1
Swiss	Q9H1C4
别名	UNC93B1;IIAE1;UNC93;UNC93B;Unc-93B1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn