

RNF125 Rabbit pAb

货号: B13909

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens,Mus musculus IF: Homo sapiens,Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	26kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,K-562
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human RNF125 (NP_060301.2).
序列	MGSVLSTDSGKSAPASATARALERRRDPELPVTSFDCAVCLEVLHQPVTRTRCGHVFCRSCIATSLKNNKWTCPCRAYLPS EGVPATDVAKRMKSEYKNCAECDT

靶点信息

研究背景	This gene encodes a novel E3 ubiquitin ligase that contains a RING finger domain in the N-terminus and three zinc-binding and one ubiquitin-interacting motif in the C-terminus. As a result of myristoylation, this protein associates with membranes and is primarily localized to intracellular membrane systems. The encoded protein may function as a positive regulator in the T-cell receptor signaling pathway.
基因ID	54941

基因名	RNF125
Swiss	Q96EQ8
别名	RNF125;TNORS;TRAC-1;TRAC1

产品验证



Western blot analysis of RNF125 expressed in Raji,K-562 using RNF125 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn