

SCNN1G Rabbit pAb

货号: B13915

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	74kDa
实测分子量	74kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,COS1,COS7,A549,NCI-H460,mouse kidney,rat kidney
细胞定位	Apical cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-250 of human SCN N1G (NP_001030.2).
序列	CNINPYKYSTVRHLLADLEQETREALKSLYGFPESRKRREAESWNSVSEGGKQPRFSHRIPLLIFDQDEKGKARDFFTGRKRK VGGSIHKASNVMHIESKQVVGFQLCSNDTSDCATYTFSSGINAIQEWYKLHYMNIMAQVPLEKKINMS

靶点信息

研究背景	Nonvoltage-gated, amiloride-sensitive, sodium channels control fluid and electrolyte transport across epit helia in many organs. These channels are heteromeric complexes consisting of 3 subunits: alpha, beta, a nd gamma. This gene encodes the gamma subunit, and mutations in this gene have been associated with Liddle syndrome.
------	---

基因ID	6340
基因名	SCNN1G
Swiss	P51170
别名	SCNN1G;BESC3;ENaCg;ENaCgamma;PHA1;SCNEG

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn