

SGMS1 Rabbit pAb

货号: B13920

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	25kDa/49kDa
实测分子量	49kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549
细胞定位	Golgi apparatus membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human SGMS1 (NP_671512.1).
序列	MKEVVYWSPKKVADWLLLENAMPEYCEPLEHFTGQDLINLTQEDFKKPPLCRVSSDNGQRLLDMIETLKMEHHLEAHKNG HANGHLNIGVDIPTPDGSFSIKIKPNGMPNGYRKEMIKIPPELERSQYPMEWGK

靶点信息

研究背景	The protein encoded by this gene is predicted to be a five-pass transmembrane protein. This gene may b e predominately expressed in brain.
基因ID	259230

基因名	SGMS1
Swiss	Q86VZ5
别名	SGMS1;MOB;MOB1;SMS1;TMEM23;hmob33

产品验证



Western blot analysis of SGMS1 expressed in A-549 using SGMS1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn