

FIP200 Rabbit pAb

货号: B13948

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	182kDa/183kDa
实测分子量	225KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U2OS,HeLa,Mouse testis,Rat testis
细胞定位	Cytoplasm,Nucleus,Peripheral membrane protein,Preautophagosomal structure membrane,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1280-1400 of human FIP 200 (NP_055596.3).
序列	TRGDSSSLVAELQEKLQEEKAKFLEQLEEQEKRKNEEMQNVRTSLIAEQQTNFNTVLTREKMRKENIINDLSDKLKSTMQQ QERQDKDLIESLSSEDRARLLEEKKKLEEEVSKLRSSSFVPS

靶点信息

研究背景	The protein encoded by this gene interacts with signaling pathways to coordinately regulate cell growth, cell proliferation, apoptosis, autophagy, and cell migration. This tumor suppressor also enhances retinobl astoma 1 gene expression in cancer cells. Alternative splicing results in multiple transcript variants encod ing distinct isoforms.
基因ID	9821

基因名	RB1CC1
Swiss	Q8TDY2
别名	RB1CC1;ATG17;CC1;FIP200;PPP1R131

产品验证



Western blot analysis of FIP200 expressed in Mouse testis,Rat testis,U2OS using FIP200 Rabbit pAb at 1 :1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn