

NDUFAB1 Rabbit pAb

货号: B13950

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rana sylvatica
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	17Da
实测分子量	12KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,NCI-H460,Mouse kidney,Mouse heart,Rat kidney,Rat heart
细胞定位	Mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-156 of human NDUFAB1 (NP_004994.1).
序列	PALVLAQVPGRVTQLCRQYSMPPLTLEGIQDRVLYVLKLYDKIDPEKLSVNSHFMKDLGLDSLQVEIIMAMEDEFGFEIPDIDAEKLMCPQEIVDYIADKKDVYE

靶点信息

研究背景	Carrier of the growing fatty acid chain in fatty acid biosynthesis (By similarity. Accessory and non-catalytic subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I, which functions in the transfer of electrons from NADH to the respiratory chain.
基因ID	4706

基因名	NDUFAB1
Swiss	O14561
别名	NDUFAB1;ACP;ACP1;FASN2A;SDAP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn