

ABCG4 Rabbit pAb

货号: B13954

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	60kDa/64kDa/65kDa/71kDa
实测分子量	72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549,U-937,Mouse brain,Rat brain
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 270-370 of human ABC G4 (NP_071452.2).
序列	KLYILSQGQCIFKGVVTNLIPYLKGLGLHCPTYHNPADFIIEVASGEYGDLPMLFRAVQNLGCAMAEKKSSPEKNEVPAPCP PCPPEVDPIESHTFATST

靶点信息

研究背景	The protein encoded by this gene is a member of the ATP-binding cassette (ABC) transporter superfamily. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). The encoded protein is a member of the White subfamily and plays an important role in cellular cholesterol homeostasis. This protein functions as either a homodimer or as a heterodimer with another ABC subfamily protein such as ABCG1.
------	--

基因ID	64137
基因名	ABCG4
Swiss	Q9H172
别名	ABCG4;WHITE2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn