

RNF6 Rabbit pAb

货号: B13960

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	16kDa/37kDa/78kDa
实测分子量	78kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	DU145,U-87MG,Mouse testis,Rat testis
细胞定位	Cell projection,Cytoplasm,Nucleus,PML body,axon
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human RNF6 (NP_898864.1).
序列	MNQSRSRSDGGSEETLPQDHNHHENERRWQQLRHREEAYYQFINELNDEDYRLMRDHNLLGTPGEITSEELQQRLDG VKEQLASQPDLRDGTNYRDSEVPRESSHEDSLLEWLNTFRRTGNATRSGQNGNQTWRAVSRTNPNNGEFRFSLEIHVN HENRGFEIHGEDYTDIPLSDSNRDHTANRQQRSTSPVARRTRSQTSVNFNGSSSNIPRTRLASRGQNPAEGSFSTLGRLR NGIGGAAGIPRANASRTNFSSTNQSGGSELRQREGQRFGAAHVWENGARSNVTVRNTNQRLPE

靶点信息

研究背景	The protein encoded by this gene contains a RING-H2 finger motif. Deletions and mutations in this gene were detected in esophageal squamous cell carcinoma (ESCC), suggesting that this protein may be a potential tumor suppressor. Studies of the mouse counterpart suggested a role of this protein in the transcription regulation that controls germinal differentiation. Multiple alternatively spliced transcript variants encoding the same protein are observed.
------	---

基因ID	6049
基因名	RNF6
Swiss	Q9Y252
别名	RNF6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn