

PHD2/EGLN1 Rabbit pAb

货号: B13963

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	36kDa/43kDa/46kDa
实测分子量	50KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,A-549,Mouse heart
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PHD2/EGLN1 (NP_071334.1).
序列	MANDSGGPGGPPSPSERDRQYCELCGKMENLLRCSRCSFFYCCKEHQRQDWKKHKLVCQGSEGALGHGVGPHQHSGPAPPAAVPPPAGAREPRKAAARR

靶点信息

研究背景	The protein encoded by this gene catalyzes the post-translational formation of 4-hydroxyproline in hypoxia-inducible factor (HIF) alpha proteins. HIF is a transcriptional complex that plays a central role in mammalian oxygen homeostasis. This protein functions as a cellular oxygen sensor, and under normal oxygen concentration, modification by prolyl hydroxylation is a key regulatory event that targets HIF subunits for proteasomal destruction via the von Hippel-Lindau ubiquitylation complex. Mutations in this gene are associated with erythrocytosis familial type 3 (ECYT3).
基因ID	54583
基因名	EGLN1
Swiss	Q9GZT9
别名	EGLN1;C1orf12;ECYT3;HALAH;HIF-PH2;HIFPH2;HPH-2;HPH2;PHD2;SM20;ZMYND6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn