

OTUD6B Rabbit pAb

货号: B13967

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Danio rerio
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:50 - 1:200
理论分子量	21kDa/33kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Raji,293T,A-431,Jurkat,Mouse brain,Mouse heart,Rat liver
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 31-323 of human OTUD 6B (NP_057107.3).
序列	MEAVLTEELDEEEQLRRHRKEKELQAKIQGMKNAVPKNDKKRRKQLTEDVAKLEKEMEQQKHREELEQLKLTTKENKIDS VAVNISNLVLENQPPRISKAQKRREKKAALEKEREERIAEAEIENLTGARHMESEKLAQILAARQLEIKQIPSDGHCMYKAIED QLKEKDCALTVVALRSQTAEYMQSHVEDFLPFLTNPNTGDMYTPEEFQKYCEDIVNTAAWGGQLELRALSHILQTPIEIIQA DSPPIIVGEEYSKKPLILVYMRHAYGLGEHYSVTRLVNIVTENC

靶点信息

研究背景	This gene encodes a member of the ovarian tumor domain (OTU)-containing subfamily of deubiquitinatin g enzymes. Deubiquitinating enzymes are primarily involved in removing ubiquitin from proteins targeted for degradation. This protein may function as a negative regulator of the cell cycle in B cells.
------	---

基因ID	51633
基因名	OTUD6B
Swiss	Q8N6M0
别名	OTUD6B;CGI-77;DUBA-5;DUBA5;IDDFSDA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn