

BRIX1 Rabbit pAb

货号: B13969

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100
理论分子量	41kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,A-549,Mouse spleen,Mouse testis,Rat spleen,Rat testis
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-353 of human BRIX1 (NP_060791.3).
序列	MAATKRKRGGFAVQAKKPKRNEIDAEPPAKRHATAEEVEEEERDRIPGPVCKGKWKNKERILIFSSRGINFRTRHLMQDL RMLMPHSKADTKMDRKDKLFVINEVCEMKNCKNCIYFEAKKKQDLYMWLSNSPHGPSAKFLVQNIHTLAELKMTGNCLK GSRPLLSFDPAFDELPHYALLKELLIQIFSTPRYHPKSQPFVDHVFTFTILDNRIWFRNFQIIIEEDAALVEIGPRFVLNLKIFQG SFGGPTLYENPHYQSPNMHRRVIRSITAACYREKQQVKDVQKLKKKEPKTLLPHDPTADVFTPAEEKPIEQWVKPEPKVD LKARKKRIYKRQRKMKQRMDSGGTK

靶点信息

研究背景	Required for biogenesis of the 60S ribosomal subunit.
基因ID	55299

基因名	BRIX1
Swiss	Q8TDN6
别名	BRIX1;BRIX;BXDC2;BRX1

产品验证



Western blot analysis of BRIX1 expressed in Mouse spleen,Rat spleen,U-87MG using BRIX1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn