

NDUFAF4 Rabbit pAb

货号: B13974

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	20kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,K-562
细胞定位	Lipid-anchor,Membrane,Mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human NDUFAF4 (NP_054884.1).
序列	MGALVIRGIRNFNLENRAEREISKMKPSVAPRHPSTNSLLREQISLYPEVKGEIARKDEKLLSFLKDVYVDSKDPVSSLQVKA AETCQEPKEFRLPKDHHFDMINIKSIPKGKISIVEALLNNHKLFPETWTAEKIMQEYQLEQKDVNSLLKYFVTFEVEIFPPED KK

靶点信息

研究背景	NADH:ubiquinone oxidoreductase (complex I) catalyzes the transfer of electrons from NADH to ubiquinone (coenzyme Q) in the first step of the mitochondrial respiratory chain, resulting in the translocation of protons across the inner mitochondrial membrane. This gene encodes a complex I assembly factor. Mutations in this gene are a cause of mitochondrial complex I deficiency.
------	---

基因ID	29078
基因名	NDUFAF4
Swiss	Q9P032
别名	NDUFAF4;C6orf66;HRPAP20;HSPC125;My013;bA22L21.1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn