

SESN2 Rabbit pAb

货号: B13982

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	54kDa
实测分子量	62kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-87MG,HepG2,K-562
细胞定位	Atg1/ULK1 kinase complex,cytoplasm,cytosol,mitochondrion,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 211-480 of human SESN 2 (NP_113647.1).
序列	VFGCGILPEGDADGSPAPQAPTPPSEQSSPPSRDPLNNSGGFESARDVEALMERMQQQLQESLLRDEGTSQEEMESRFEL EKSESLLVTPSADILEPSHPDMLCFVEDPTFGYEDFTRRGAQAPPTFRAQDYTWEDHGYSLIQRLYPEGGQLLDEKFQAA YSLTYNTIAMHSGVDTSVLRRAIWNYIHCVFGIRYDDYDGEVNQLLERNLKVYIKTVACYPEKTTRRMYNLFWRHFRHSEK VHVNLLLLEARMQAALLYALRAITRYMT

靶点信息

研究背景	This gene encodes a member of the sestrin family of PA26-related proteins. The encoded protein may function in the regulation of cell growth and survival. This protein may be involved in cellular response to different stress conditions.
基因ID	83667

基因名	SESN2
Swiss	P58004
别名	SESN2;HI95;SES2;SEST2;sestrin-2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn