

IGFBP5 Rabbit pAb

货号: B14003

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rattus norvegicus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	30kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 21-272 of human IGFBP 5 (NP_000590.1).
序列	LGSFVHCEPCDEKALSMCPPSPLGCELVKEPGCGCCMTCALAEGQSCGVYTERCAQGLRCLPRQDEEKPLHALLHGRGVCLNEKSYREQVKIERDSREHEEPTTSEMAEETYSPIFRPKHTRISELKAEAVKKDRRKLTQSKFVGGAENTAHPRIISAPEMRQESEQGPCRRRHMEASLQELKASPRMVPRAVYLPNCDRKGIFYKRKQCKPSRGRKRGICWCVDKYGMKLPGMEYVDGDFQCHTFDSSNVE

靶点信息

研究背景	IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.
------	---

基因ID	3488
基因名	IGFBP5
Swiss	P24593
别名	IGFBP5;IBP5

产品验证



Western blot analysis of IGFBP5 expressed in human lung cancer, rat testis using IGFBP5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn