

DUSP12 Rabbit pAb

货号: **B14048**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	37kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,A-549
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 130-270 of human DUSP12 (NP_009171.1).
序列	LMKTDQLPFEKAYEKLQILKPEAKMNEGFEWQLKLYQAMGYEVDTS SAIYKQYRLQKVTEKYPELQNL PQELFAVDPTTVS QGLKDEVLYKCRKRRSLFRSS SILDHREGSGPIAFAHKRMT PSSMLTTGRQAQCTSYFI

靶点信息

研究背景	The protein encoded by this gene is a member of the dual specificity protein phosphatase subfamily. These phosphatases inactivate their target kinases by dephosphorylating both the phosphoserine/threonine and phosphotyrosine residues. They negatively regulate members of the mitogen-activated protein (MAP) kinase superfamily (MAPK/ERK, SAPK/JNK, p38), which is associated with cellular proliferation and differentiation. Different members of the family of dual specificity phosphatases show distinct substrate specificities for various MAP kinases, different tissue distribution and subcellular localization, and different modes of inducibility of their expression by extracellular stimuli. This gene product is the human ortholog of the <i>Saccharomyces cerevisiae</i> YVH1 protein tyrosine phosphatase. It is localized predominantly in the nucleus, and is novel in that it contains, and is regulated by a zinc finger domain.
基因ID	11266
基因名	DUSP12
Swiss	Q9UNI6
别名	DUSP12;DUSP1;YVH1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn