

RBM14 Rabbit pAb

货号: B14056

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:100 - 1:500
理论分子量	12kDa/13kDa/17kDa/37kDa/69kDa
实测分子量	75KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,Jurkat,Rat testis
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-220 of human RBM 14 (NP_006319.1).
序列	VVKDYAFVHMEKEADAKAAIAQLNGKEVKGKRINVELSTKGQKKGPGGLAVQSGDKTKKPGAGDTAFPGTGGFSATFDYQ QAFGNSTGGFDGQARQPTPPFFGRDRSPLRRS

靶点信息

研究背景	This gene encodes a ribonucleoprotein that functions as a general nuclear coactivator, and an RNA splicing modulator. This protein contains two RNA recognition motifs (RRM) at the N-terminus, and multiple hexapeptide repeat domain at the C-terminus that interacts with thyroid hormone receptor-binding protein (TRBP), and is required for transcription activation. Alternatively spliced transcript variants encoding different isoforms (with opposing effects on transcription) have been described for this gene.
------	---

基因ID	10432
基因名	
Swiss	Q96PK6
别名	RBM14;COAA;PSP2;SIP;SYTIP1;TMEM137

产品验证



Western blot analysis of RBM14 expressed in Rat testis, HeLa, 293T using RBM14 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn