

HAL Rabbit pAb

货号: B14065

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Helicoverpa armigera
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	49kDa/64kDa/72kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HAL (N P_002099.1).
序列	MPRYTVHVRGEWLAVPCQDAQLTVGWLGREAVRRYIKNKPDNGGFTSVDDAHFLVRRCKGLGLLDNEDRLEVALENNE FVEVVIEGDAMSPDFIPSQPEGVYLYSKYREPEKYIELDGDRLTTEDLVNLGKGRYKIKLTPTAEKRVQKSREVIDSIIKEKTV VYGITTGFGKFARTVIPINKLQELQVNLVRSRSSGVGK

靶点信息

研究背景	Histidine ammonia-lyase is a cytosolic enzyme catalyzing the first reaction in histidine catabolism, the no noxidative deamination of L-histidine to trans-urocanic acid. Histidine ammonia-lyase defects cause histi dinemia which is characterized by increased histidine and histamine and decreased urocanic acid in body fluids. Several transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	3034
基因名	HAL
Swiss	P42357
别名	HAL;HIS;HSTD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn