

ARR3 Rabbit pAb

货号: B14068

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Cynoglossus semilaevis
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa/42kDa
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse eye,Rat eye
细胞定位	cytoplasm,photoreceptor inner segment,photoreceptor outer segment
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ARR3 (NP_004303.2).
序列	MSKVFKKTSSNGKLSIYLGKRDFVDHVDTVEPIDGVVLVDPEYLKCRKLFVMLTCAFRYGRDDLEVIGLTFRKDLYVQTLQ VVPAESSSPQGPLTVLQERLLHKLGDNAYPFTLQMVTNLPCSVTLQPGPEDAGKPCGIDFEVKSFCAENPEETVSKRDYVR LVVRKVQFAPPEAGPGPSAQTIIRFLLSAQPLQLQAWMDREVHYHGEPISVNVSINNCTNKVIKKIKISVDQITDVVLYSLDK YTKTVFIQEFTETVA

靶点信息

研究背景	The protein encoded by this gene is a non-visual arrestin which binds to agonist-activated, phosphorylate d G protein-coupled receptors. This binding uncouples the receptor from the heterotrimeric G protein, res ulting in termination of the G protein-coupled receptor signaling. The encoded protein also is a part of the centrosome, interacting with gamma-tubulin to help regulate proper centrosome function.
------	--

基因ID	407
基因名	ARR3
Swiss	P36575
别名	ARR3;ARRX;cArr

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn