

RAB14 Rabbit pAb

货号: B14086

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	23kDa
实测分子量	24KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Neuro-2a
细胞定位	Cytoplasmic side,Cytoplasmic vesicle,Early endosome membrane,Golgi apparatus,Golgi apparatus membrane,Lipid-anchor,Recycling endosome,phagosome,phagosome membrane,trans-Golgi network membrane
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human RAB14 (NP_057406.2).
序列	TRIIEVSGQIKIKLIWDTAGQERFRAVTRSYRGAAGALMVYDITRRSTYNHLSSWLTDARNLTNPNTVIILIGNKADLEAQRDVTYEEAKQFAEENGLLF

靶点信息

研究背景	RAB14 belongs to the large RAB family of low molecular mass GTPases that are involved in intracellular membrane trafficking. These proteins act as molecular switches that flip between an inactive GDP-bound state and an active GTP-bound state in which they recruit downstream effector proteins onto membranes (J unutula et al., 2004 [PubMed 15004230]).
------	---

基因ID	51552
基因名	RAB14
Swiss	P61106
别名	RAB14;FBP;RAB-14

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn