

AIFM2 Rabbit pAb

货号: B14128

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	36kDa/40kDa
实测分子量	41KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse liver,Mouse testis,Rat liver,Rat testis
细胞定位	Cytoplasm,Lipid-anchor,Membrane,Mitochondrion outer membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 43-371 of human AIFM2 /AMID/AMID (NP_116186.1).
序列	KDSFHHNVAALRASVETGFAKKTFFISYSVTFKDNFRQGLVVGIDLKNQMVLLQGGEALPFSHLILATGSTGPFPGKFNEVS SQQAAIQAYEDMVRQVQRSRFIVVVGGSAGVEMAAEIKTEYPEKEVTLIHSQVALADKELLPSVRQEVKEILLRKGVQLLL SERVSNLEELPLNEYREYIKVQTDKGTEVATNLVILCTGIKINSSAYRKAFESRLASSGALRVNEHLQVEGHSNVYAIGDCAD VRTPKMAYLAGLHANIAVANIVNSVKQRPLQAYKPGALTFLLSMGRNDGVGQISGFYVGRMLMVRLTCSRDLFVSTSWKTM RQS

靶点信息

研究背景	This gene encodes a flavoprotein oxidoreductase that binds single stranded DNA and is thought to contri bute to apoptosis in the presence of bacterial and viral DNA. The expression of this gene is also found to be induced by tumor suppressor protein p53 in colon cancer cells.
------	--

基因ID	84883
基因名	AIFM2
Swiss	Q9BRQ8
别名	AIFM2;AMID;PRG3; FSP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn