

CLEC2D Rabbit pAb

货号: B14142

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	10kDa/14kDa/17kDa/21kDa/22kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,HL-60,LO2,U-87MG,Mouse heart
细胞定位	Cell membrane,Cell surface,Endoplasmic reticulum,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-154 of human CLEC2 D (NP_037401.1).
序列	RANCHQEPSVCLQAACPESWIGFQRKCFYFSDDTKNWTSSQRFCDSDADLAQVESFQELNFLRLRYKGPSDHWIGLSR EQGQPWKWINGTEWTRQ

靶点信息

研究背景	This gene encodes a member of the natural killer cell receptor C-type lectin family. The encoded protein i nhibits osteoclast formation and contains a transmembrane domain near the N-terminus as well as the C-type lectin-like extracellular domain. Several alternatively spliced transcript variants have been identified for this gene.
基因ID	29121

基因名	CLEC2D
Swiss	Q9UHP7
别名	CLEC2D;CLAX;LLT1;OCIL

产品验证



Western blot analysis of CLEC2D expressed in Mouse heart,Raji,HL-60 using CLEC2D Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn