

# LIMA1 Rabbit pAb

货号: B14162

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                          |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  | <b>WB:</b> Mus musculus                                                                              |
| 应用    | <a href="#">WB</a> <a href="#">IHC</a>                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:100 - 1:200                                               |
| 理论分子量 | 51kDa/67kDa/85kDa                                                                                    |
| 实测分子量 | 100kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | A-549                                                                                                |
| 细胞定位  | Cell junction,Cytoplasm,cytoskeleton,focal adhesion                                                  |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                                                                                                |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 560-759 of human LIMA 1 (NP_057441.1).                                                                                           |
| 序列   | DEISKPEVPEDVDLDLKKLRRSSSLKERSRPFTVAASFQSTSVKSPKTVSPPIRKGWSMSEQSEESVGGRVAERKQVENA<br>KASKKNGNVGKTTWQNKESKGETGKRSGEGHSLEMENENLVENGADSDDDNSFLKQQSPQEPKSLNWSSFVDNTFA<br>EEFTTQNKSQDVELWEGEVVKELSVEEQIKRNRYYDEDEDEE |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                        |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a cytoskeleton-associated protein that inhibits actin filament depolymerization and cross-links filaments in bundles. It is downregulated in some cancer cell lines. Alternatively spliced transcript variants encoding different isoforms have been described for this gene, and expression of some of the variants may be independently regulated. |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                    |
|-------|--------------------|
| 基因ID  | 51474              |
| 基因名   | LIMA1              |
| Swiss | Q9UHB6             |
| 别名    | LIMA1;EPLIN;SREBP3 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)