

GAS2 Rabbit pAb

货号: B14163

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Oreochroms mossambcus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	15kDa/34kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,Mouse liver,Mouse kidney
细胞定位	Cytoplasm,Membrane,Peripheral membrane protein,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human GAS2 (NP_005247.1).
序列	MCTALSPKVRSGPGLSDMHQYSQWLASRHEANLLPMKEDLALWLTNLLGKEITAETFMEKLDNGALLCQLAETMQEKFK ESMDANKPTKNLPLKKIPCKT

靶点信息

研究背景	The protein encoded by this gene is a caspase-3 substrate that plays a role in regulating microfilament an d cell shape changes during apoptosis. It can also modulate cell susceptibility to p53-dependent apoptosi s by inhibiting calpain activity. Alternative splicing results in multiple transcript variants.
基因ID	2620

基因名	GAS2
Swiss	O43903
别名	GAS2;GAS-2

产品验证



Western blot analysis of GAS2 expressed in SW620,Mouse liver,Mouse kidney using GAS2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn