

MYT1 Rabbit pAb

货号: B14191

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	122kDa/125kDa
实测分子量	122kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,U-251MG,Mouse liver,Mouse kidney,Rat kidney,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 650-820 of human MYT 1 (NP_004526.1).
序列	KPQDLPSKSVDIEVDENGLDLSMHKHKRENAPSSSSCSSSPGVKSPDASQRHSSTSAPSSSMTSPQSSQASRQDEWDRPLDYTKPSRLREEEP EESEPAAHSFASSEADDQEVSEENFEERKYPGEVTLTNFKLKFLSKDIKKELLCPTPGCDGSGHITGNYASHRS

靶点信息

研究背景	The protein encoded by this gene is a member of a family of neural specific, zinc finger-containing DNA-binding proteins. The protein binds to the promoter regions of proteolipid proteins of the central nervous system and plays a role in the developing nervous system.
基因ID	4661

基因名	MYT1
Swiss	Q01538
别名	MYT1;C20orf36;MTF1;MYTI;NZF2;PLPB1;ZC2H2C1;ZC2HC4A

产品验证



Western blot analysis of MYT1 expressed in Mouse liver,Rat kidney,U-87MG using MYT1 Rabbit pAb at 1 :1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn