

# TRIM31 Rabbit pAb

货号: B14202

## 产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                            |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  | <b>IHC:</b> Mus musculus                                                                               |
| 应用    | <a href="#">WB</a> <a href="#">IF/ICC</a>                                                              |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                               |
| 理论分子量 | 30kDa/48kDa                                                                                            |
| 实测分子量 | 48kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | HT-29,Mouse large intestine,Mouse small intestine                                                      |
| 细胞定位  | Cytoplasm,Mitochondrion                                                                                |
| 纯化    | Affinity purification                                                                                  |

## 抗原信息


|      |                                                                                                                                                                                                                                                        |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 186-425 of human TRIM 31 (NP_008959.3).                                                                                                                                  |
| 序列   | ELLHQVLEEEKNFLLSRIYWLGHEGTEAGKHYVASTEPQLNDLKKLVDSLKTKQNMPPRQLLEDIKVVLCRSEEFQFLNPTP<br>VPLELEKKLSEAKSRHDSITGSLKKFKDQLQADRKKDENRFFKSMNKNDMKSWGLLQKNNHKMNKTSEPGSSSAGGRT<br>TSGPPNHHSAPSHSLFRASSAGKVTFPVCLASYDEISGQGASSQDTKTFDVALUEELHAALSEWLTAIRAWFCEVPSS |

## 靶点信息

|      |                                                                                                                                                                                                                                                  |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a protein that functions as an E3 ubiquitin-protein ligase. This gene shows altered expression in certain tumors and may be a negative regulator of cell growth. Alternative splicing results in multiple transcript variants. |
| 基因ID | 11074                                                                                                                                                                                                                                            |

|       |                              |
|-------|------------------------------|
| 基因名   | TRIM31                       |
| Swiss | Q9BZY9                       |
| 别名    | TRIM31;C6orf13;HCG1;HCGI;RNF |

产品验证



Western blot analysis of TRIM31 expressed in HT-29,Mouse large intestine,Mouse small intestine using TRIM31 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)