

Arp3 Rabbit pAb

货号: B14203

产品信息

反应	Human,Mouse,Rat,Monkey
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:10 - 1:100
理论分子量	47kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,HL-60,HeLa,Jurkat,COS-1,COS-7,Mouse spleen
细胞定位	Cell projection,Cytoplasm,cytoskeleton
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-418 of human Arp3 (N P_005712.1).
序列	MAGRLPACVVDCGTGYTKLGYAGNTEPQFIIPSCIAIKESAKVGDDAQRRVMKGVDDLDFFIGDEAIEKPTYATKWPIRHGI VEDWDLMERFMEQVIFKYLRAEPEDHYFLLTEPLNTPENREYTAEIFESFNVPGLYIAVQAVLALAASWTSRQVGERTLT GTVIDSGDGVTHVIPVAEGYVIGSCIKHIPIAGRDITYFIQQLLRDREVGIPPEQSLETAKAVKERYSYVCPDLVKEFNKYDTD GSKWIKQYTGINAISKKEFSIDVGYERFLGPEIFFHPEFANPDFTQPISEVVDEVIQNCPIDVRRPLYKNIVLSGGSTMFRDFG RRLQRDLKRTVDARLKLSEELSGGRLKPKPIDVQVITHHMQRVAVWFGGSMLASTPEFYQVCHTKKDYEEIGPSICRHNPV FGVMS

靶点信息

研究背景	The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Three transcript variants encoding two different isoforms have been found for this gene.
基因ID	10096
基因名	ACTR3
Swiss	P61158
别名	ACTR3;ARP3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn