

SEC11A Rabbit pAb

货号: B14210

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Homo sapiens IHC: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	17kDa/18kDa/20kDa/21kDa
实测分子量	21kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SGC-7901,HT-29,Raji,HepG2,A-549,Mouse kidney,Mouse liver,Mouse lung,Rat brain
细胞定位	Endoplasmic reticulum membrane,Microsome membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 37-145 of human SEC11A (NP_055115.1).
序列	KGLMVTGSESPIVVVLSGSMEPAFHRGDLLFTNRPVEDPIRVGEIVVFRIEGREIPIVHRVLKIHEKQNGHIKFLTKGDNNAVDDRGLYKQGQHWLEKKDVVGRARGF

靶点信息

研究背景	This gene encodes a member of the peptidase S26B family. The encoded protein is an 18kDa subunit of the signal peptidase complex and has been linked to cell migration and invasion, gastric cancer and lymph node metastasis. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 8.
------	--

基因ID	23478
基因名	SEC11A
Swiss	P67812
别名	SEC11A;1810012E07Rik;SEC11L1;SPC18;SPCS4A;sid2895

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn