

NOXO1 Rabbit pAb

货号: B14217

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	40kDa/41kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse small intestine,Mouse stomach
细胞定位	Cell membrane,Cytoplasmic side,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 121-370 of human NOX O1 (NP_653204.1).
序列	QPLDLEPALPPGSRVILPTPEEQPLSRAAGRLSIHSLEAQLRCLQPFCQDTRDRPFQAQAQESLDVLLRHPSGWWLVEN EDRQTAWFPAPYLEEAAPGQGREGGPSLGSSGPQFCASRAYESSRADELSVPAGARVRVLETSDRGWWLCRYGDRAGL LPAVLLRPEGLGALLSGTGFRGGDDPAGEARGFPEPSQATAPPPTVPTRPSPGAIQSRCCTVTRRALERRPRRQGRPRGCV DSVPHPTTEQ

靶点信息

研究背景	This gene encodes an NADPH oxidase (NOX) organizer, which positively regulates NOX1 and NOX3. The p rotein contains a PX domain and two SH3 domains. Alternatively spliced transcript variants encoding mult iple isoforms have been observed for this gene.
基因ID	124056

基因名	NOXO1
Swiss	Q8NFA2
别名	NOXO1;P41NOX;P41NOXA;P41NOXB;P41NOXC;SH3PXD5;SNX28

产品验证



Western blot analysis of NOXO1 expressed in Mouse small intestine, Mouse stomach using NOXO1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn