

RNF138 Rabbit pAb

货号: B14227

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Sus scrofa
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	24kDa/28kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-245 of mouse RNF138 (Q9CQE0).
序列	SVTRRERACPERALDLENIMRRFSGSCRCCKIKFYRMRHHYKSCCKYQDEYGVSSVIPNFKISQDSVRSSNRSETSASDNTETYQEDTSSSGHPTFKCPLCQESNFTQRLLDHCNSNHLFQIVPVTCPICVSLPWGDP SQITRNFVSHLNQRHQFDYGEFVNLQLDEETQYQTAVEESFQVNM

靶点信息

研究背景	The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. Alternatively spliced transcript variants encoding distinct isoforms have been observed.
基因ID	56515

基因名	
Swiss	
别名	RNF138;HSD-4;NARF;STRIN;hNARF

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn