

CDK13 Rabbit pAb

货号: B14234

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	158kDa/164kDa
实测分子量	164kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG
细胞定位	Nucleus speckle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 620-900 of human CDK 13 (NP_003709.3).
序列	KEADSLRGNISVKAVKKEVEKKLRCLLADLPLPELPGDDLSKSPEEKKTATQLHSKRRPKICGPRYGETKEKDIDWGKRC VDKFDIIGIIGEGTYGQVYKARDKDTGEMVALKKVRLDNEKEGFPITAIREIKILRQLTHQSIINMKEIVTDKEDALDFKKDKG AFYLVFEYMDHDLMLLESGLVHFNENHIKSFMRQLMEGLDYCHKKNFLHRDIKCSNILLNNRGQIKLADFGFLARLYSSEE SRPYTNKVITLWYRPPPELLLGEERYTPAIDVWSC

靶点信息

研究背景	The protein encoded by this gene is a member of the cyclin-dependent serine/threonine protein kinase family. Members of this family are well known for their essential roles as master switches in cell cycle control. The exact function of this protein has not yet been determined, but it may play a role in mRNA processing and may be involved in regulation of hematopoiesis. Alternatively spliced transcript variants have been described.
------	---

基因ID	8621
基因名	CDK13
Swiss	Q14004
别名	CDK13;CDC2L;CDC2L5;CHDFIDD;CHED;hCDK13

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn