

PLC beta 3 (PLCB3) Rabbit pAb

货号: B14241

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Rattus norvegicus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	131kDa/138kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR8,Mouse lung,Rat lung
细胞定位	Membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1010-1234 of human PLC beta 3 (PLCB3) (NP_001303243.1).
序列	GALGGAADVEDTKEGEDEAKRYQEFQNRQVQSLLLELREAQVDAEAQRRLEHLRQALQRLREVLDANTTQFKRLKEMNE REKKELQKILDRKRHNSISEAKMRDKHKKEAELTEINRRHITESVNSIRRLEEAQQRHDLVAGQQQVLQQLAEEEPKLLA QLAQECQEQRARLPQEIRRSLLGEMPEGLGDGPLVACASNGHAPGSSGHLSGADSESQEENTQL

靶点信息

研究背景	This gene encodes a member of the phosphoinositide phospholipase C beta enzyme family that catalyze the production of the secondary messengers diacylglycerol and inositol 1,4,5-triphosphate from phosphatidylinositol in G-protein-linked receptor-mediated signal transduction. Alternative splicing results in multiple transcript variants.
------	--

基因ID	5331
基因名	PLCB3
Swiss	Q01970
别名	PLCB3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn