

MFAP2 Rabbit pAb

货号: B14243

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	17kDa/20kDa
实测分子量	31kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,OVCAR3,MCF7,Mouse lung,Mouse brain,Rat lung
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 18-183 of human MFAP2 (NP_059453.1).
序列	QGQYDLDPFPDPDHVQYTHYSDQIDNPDDYDYQEVTPRPSEEQFQFQSQQVQQEVIPAPTPEPGNAELEPTTEPGPLDC REEQYPCTRLYSIHRPCKQCLNEVCFYSLRRVYVINKEICVRTVCAHEELLRADLCRDKFSKCGVMASGLCQSVAASCARS CGSC

靶点信息

研究背景	Microfibrillar-associated protein 2 is a major antigen of elastin-associated microfibrils and a candidate for involvement in the etiology of inherited connective tissue diseases. Four transcript variants encoding two different isoforms have been found for this gene.
基因ID	4237

基因名	MFAP2
Swiss	P55001
别名	MFAP2;MAGP;MAGP-1;MAGP1

产品验证



Western blot analysis of MFAP2 expressed in Mouse lung,Rat lung,HeLa using MFAP2 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn