

GOLGA4 Rabbit pAb

货号: B14244

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:1000 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	260kDa/261kDa/262kDa
实测分子量	280kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,SW480,A-549,MCF7,Mouse heart,Mouse brain
细胞定位	Cytoplasm,Golgi apparatus,Golgi apparatus membrane,Peripheral membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human GOLGA 4 (NP_002069.2).
序列	MFKKLKQKISEEQQLQQALAPAQASSNSSTPTRMRSRTSSFTEQLDEGTPNRESGDTQSFAQKLQLRVPSVESLFRSPI KESLFRSSSKESLVRTSSRESLNRLDLSSTASFDPSPDMDSEAEDLVGNSDSLNKEQLIQLRRMERSL

靶点信息

研究背景	The Golgi apparatus, which participates in glycosylation and transport of proteins and lipids in the secretory pathway, consists of a series of stacked cisternae (flattened membrane sacs). Interactions between the Golgi and microtubules are thought to be important for the reorganization of the Golgi after it fragments during mitosis. This gene encodes one of the golgins, a family of proteins localized to the Golgi. This protein has been postulated to play a role in Rab6-regulated membrane-tethering events in the Golgi apparatus. Alternatively spliced transcript variants encoding different isoforms have been identified in this gene.
基因ID	2803
基因名	GOLGA4
Swiss	Q13439
别名	GOLGA4;CRPF46;GCP2;GOLG;MU-RMS-40.18;p230;golgin A4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn