

LRP4 Rabbit pAb

货号: B14250

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	212kDa
实测分子量	200kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,HepG2,Mouse lung,Mouse liver
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1500-1725 of human LR P4 (NP_002325.2).
序列	NLDGSEKVLINTDLGWPNGLTLDYDTRRIYWVDAHLDRIESADLNGKLRQVLVSHVSHPFALTQQDRWIYWTDWQTKS IQRVDKYSGRNKETVLANVEGLMDIIVVSPQRQTGTNACGVNNGGCTHLCFARASDFVCACPDEPDSRPCSLVPGLVPPA PRATGMSEKSPVLPNTPTTLYSSTTRTRTSLEEVEGRCSERDARLGLCARSNDAPPAAPGEGLHIS

靶点信息

研究背景	This gene encodes a member of the low-density lipoprotein receptor-related protein family. The encoded protein may be a regulator of Wnt signaling. Mutations in this gene are associated with Cenani-Lenz synd rome.
基因ID	4038

基因名	LRP4
Swiss	O75096
别名	LRP4;CLSS;CMS17;LRP-4;LRP10;MEGF7;SOST2

产品验证



Western blot analysis of LRP4 expressed in Mouse lung,SKOV3,HepG2 using LRP4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn