

IFNLR1 Rabbit pAb

货号: B14257

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:1000 - 1:4000
理论分子量	23kDa/27kDa/54kDa/57kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Jurkat,MCF7,Mouse lung,Mouse heart,Mouse small intestine,Rat lung,Rat heart
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 21-228 of human IFNLR 1 (NP_734464.1).
序列	RPRLAPPQNVTLLSQNFSVYLTWLPGLGNPQDVTFVAYQSSPTRRRWREVEECAGTKELLCSMMCLKKQDLYNKFKGR VRTVSPSSKSPWVESEYLDYLFVEVAPPVLVLTQTEEILSANATYQLPPCMPPLDLKYEVAFWKEGAGNKTLFPVTPHGQP VQITLQPAASEHHCLSARTIYTFSPKYSKFSKPTCFLLEVPEANWA

靶点信息

研究背景	The protein encoded by this gene belongs to the class II cytokine receptor family. This protein forms a receptor complex with interleukine 10 receptor, beta (IL10RB). The receptor complex has been shown to interact with three closely related cytokines, including interleukin 28A (IL28A), interleukin 28B (IL28B), and interleukin 29 (IL29). The expression of all three cytokines can be induced by viral infection. The cells over expressing this protein have been found to have enhanced responses to IL28A and IL29, but decreased response to IL28B. Three alternatively spliced transcript variants encoding distinct isoforms have been reported.
基因ID	163702
基因名	IFNLR1
Swiss	Q8IU57
别名	IFNLR1;CRF2/12;IFNLR;IL-28R1;IL28RA;LICR2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn